

SOLAX Genie Automatic Folding Mobility Scooter



Features

- Automatic folding/unfolding using a remote control
- Able to automatic fold/unfold using the fold/unfold switch at the back of the scooter in the event the remote control runs out of battery
- Can be manually folded/unfolded easily in the event both the remote control and scooter run out of battery
- Fits in any trunk
- Light weight, easy carry
- 24 V/10 AH lithium battery allows travel up to 15 km with option of spare battery to travel up to 30 km
- Ideal and approved for airline travel

DISTRIBUTED BY:

Independent Living Specialists

<https://ilsau.com.au>



Automatic Folding/Unfolding

Folds and unfolds automatically using a remote control or the fold/unfold switch at the back of the scooter. Can also be manually folded/unfolded easily with no tools required



Compact Size & Convenient to Travel

Light weight, easy to carry and fits in any trunk, and can be pulled around like a luggage

SPECIFICATIONS

SOLAX GENIE AUTOMATIC FOLDING MOBILITY SCOOTER

ITEM #	SOLAX-Genie
Unfolded Dimensions (L X W X H)	930 mm (L) X 450 mm (W) X 857 mm (H)
Folded Dimensions (L X W X H)	365 mm (L) X 450 mm (W) X 580 mm (H)
Seat Height	560 mm
Tire	6 inches for front wheel (solid wheel), 7 inches for rear wheel (solid wheel)
Maximum Speed	6 kmh/3.75 mph
Safe Gradient	Up to 15°
Maximum Gradient	8°
Range	15 km/9.32 miles
Turning Circle	1200 mm
Total Weight (Battery Included)	25 kg (55lbs)
Brake	Intelligent regenerative electromagnetic brake
Drive System	Rear wheel drive
Handlebar	3 positions
Maximum Capacity	120 kg
Ground Clearance	115 mm (4.5 inches)
Battery	24 V/10 AH (Lithium Battery)
Motor	120 W
Battery Charger	DC 24 V/2 A off-board (automatic protection charger)
Colors	Blue
Warranty	1 Year on Frame, Electronics and Motor 6 Months on Battery

DISTRIBUTED BY:

Independent Living Specialists

Unit 1 12 Mars Road, Lane Cove, 2066

<https://ilsau.com.au>